



PAT Conference 2008

Exercising Real Control Over API & Excipient Isolation

The Power of Sonocrystallization!

30 October 2008, Heidelberg, Germany

Christian Jones, Business Development Manager, Prosonix Ltd.

© European Compliance Academy (ECA)



Agenda

- Company Overview
- The Critical Importance of Controlled Crystallization in Manufacture
- Sonocrystallization and Ultrasonic Particle Engineering
 - The Key to Controlled Crystallization of Bulk API and Excipient
 - #1 Theory and Examples
 - #2 Commercial scale solutions
- Advanced Ultrasonic Particle Engineering of Difficult to Manufacture products, such as those for Inhalation

© European Compliance Academy (ECA)





Prosonix Ltd




- Based in Oxford, UK
- Focused on patented ultrasonic particle engineering technology using proprietary sonocrystallization techniques to make better medicines
- Income stream from partner funded collaborations, licensing, and product supply
- Development collaborations with 8 of top 10 pharmaceutical companies
- *NEW* recent Milestone success with Pfizer

© European Compliance Academy (ECA)



Prosonix Vision

"Making Better Medicines More Efficiently"



© European Compliance Academy (ECA)





The Critical Importance of Controlled Crystallization in Manufacture

Pharmaceutical Manufacturing & QbD



- **Selection Criteria for NCE's are:**
 - Clinical Efficacy
 - Bioavailability
 - Stability
 - Processability
- **Manufacturing and ultimate product success is dictated by:**
 - Product physical form
 - *polymorphism, crystal habit etc.*
 - Process scale-up and manufacturing problems
 - *raw material variations*
 - *synthetic complexity*
- Improved understanding of manufacture and design space could also lead to better asset utilisation, quality and lower failure rates
- Sonocrystallization approaches can transform productivity at the same time as improving flexibility and manufacturing compliance

Importance of Controlled Crystallization

- Crystallization is a ubiquitous and critical manufacturing unit operation
- Almost every chemical process that produces solid form involves at least one crystallization step, either for intermediate separation, final product purification, or for the removal of key impurities
- Crystallization processes are poorly understood and are difficult to control
- Control of the nucleation difficult but is the key to process control
- Process robustness governs process productivity and economics
- Physical form dictates drug product quality and effectiveness

Importance of Materials Science


A lack of understanding in early development stores up future problems



- Classic series based approach to isolation and development brings problems

Importance of Controlled Crystallization

But by fixing things early by design enhances productivity and control




- Controlled crystallisation eliminates downstream problems
- Sonocrystallization potentially “best available technology” for API manufacture

QbD Manufacture: A Cultural Change

Acknowledgement to Prof. Kevin Roberts, University of Leeds


Historic Approach


“Process Down”



Future QbD

“Molecule Up”





Sonocrystallization



The Key to Controlled Crystallization of Bulk API's and Excipients

#1 Theory and Examples




Key Benefits of Sonocrystallization

- Control particle size, shape, crystallinity, polymorphism
- Improve batch consistency, filtration, isolation and drying
- Improve formulation consistency, stability and performance
- Enhance dissolution of poorly soluble drugs
- Replace problem physical seeding
- Increase cGMP compliance



Comparison of Acoustical Energy to Other Sources



- Utilisation of ultrasound is a process intensification method in comparison to other energy applications

- Enables delivery of large quantities of energy in short timeframes on microscales, but for controlled macro process effect


Suslick, K. S.; Didenko, Y.; Fang, M. M.; Hyeon, T.; Kolbeck, K. J.; McNamara, W. B. III; Mdeleeni, M. M.; Wong, M. "Acoustic Cavitation and Its Chemical Consequences" Phil. Trans. Roy. Soc. London A, 1999, 357, 335-353.

© European Compliance Academy (ECA)

Sonocrystallization

Consistent crystallization for API and Excipients

- Classically **Nucleation** is random, and resultant crystallization processes are uncontrolled, leading to poorly performing API, and drug formulations
- Molecules of product assemble in clusters. The clusters will progressively increase in size to become viable crystals



- **Controlled** nucleation is fundamental to crystallization control
- Power **ultrasound** via cavitation allows controlled nucleation, i.e. *Sonocrystallization*
- Increasingly recognised and used for manufacturing improvement by *Merck* (Aaron Moment, Jan 2008, Leeds RSC SonoChemistry Conference), *AZ* (Literature), *Pfizer* (Material Science Presentation, Ivan Marziano APS meeting 2007)

© European Compliance Academy (ECA)

ECA **Control of Crystal Size**
 EUROPEAN COMPLIANCE ACADEMY General Rules on the Effects of Cavitation on a Saturated Solution


ULTRASONIC TREATMENT OF SUPERSATURATED SOLUTION

1. Continuous insonation produces many nuclei resulting in small crystals
2. Initial insonation produces finite nuclei which can be grown into large crystals
3. Pulsed insonation gives tailored crystal size

ULTRASONIC TREATMENT BEFORE & AFTER CRYSTALLIZATION

4. Continuous insonation throughout supersaturation produces many nuclei resulting in small crystals. Application of ultrasound thereafter can condition the crystals produced

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®



ECA **Polymorph Control with Ultrasound #1**
 EUROPEAN COMPLIANCE ACADEMY L-Glutamic Acid

- L-glutamic acid has two polymorphic forms: α & β
- Metastable α -form: produced under kinetic control
- The transformation of form α to β is solution mediated
- The metastable a form is difficult to obtain
- Use power ultrasound to reproducibly prepare the α or β form
- Sononucleation at different supersaturation
- **PAT Tool: Online Raman**


© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®

ECA **Polymorph Control with Ultrasound #2**
 EUROPEAN COMPLIANCE ACADEMY Advanced small molecule API

Key results: reduced median crystal size, tighter CSD, polymorph control
PAT Tools: Online Raman and Lasentec

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®

ECA APIs - PRT Improves flowability
 EUROPEAN COMPLIANCE ACADEMY Diltiazem Hypertension Product




Raw powder → **Rounded powder**

- Improved uniform packing density of powders
- Enhanced flow of powders
- Reduced electro-static charges
- Manufacturing by direct compression without granulation
- Higher filler loading in composite pastes
- Bypass patents on particle size/shape
- Significantly enhanced flowability
- Improved stability on storage (no agglomeration)

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®

ECA Excipients – PRT applied to Lactose
 EUROPEAN COMPLIANCE ACADEMY




Raw powder → **Rounded powder**

PAT Tool: Online Lasentec

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®

Sonocrystallization, Sonomilling, Particle rounding



Specific angiotensin II type 1 antagonist - Hypertension




Cooling crystallisation




Sonomilling and rounding



Sonocrystallisation



Sonocrystallization,
Sonomilling and rounding



- Agglomeration and solvent inclusion is avoided
- Circumvent micronisation

© European Compliance Academy (ECA)

Sonocrystallization The Key to Controlled Crystallization of Bulk API's and Excipients #2 Commercial Scale Solutions

© European Compliance Academy (ECA)

ECA
EUROPEAN COMPLIANCE ACADEMY

Probe Problems:

Cavitation surface erosion hitherto prohibited commercial use

Liquid jet penetrates bubble during asymmetric collapse

Damage to a solid caused by jet impact and emission of shock waves as a result of repetitive bubble implosions

Vessel wall
Probe tip
Solid surface
Jet impact

Bubble collapse
Bubble implosion
Liquid in-flow

Probes can not be used for commercial scale production

© European Compliance Academy (ECA)

Prosonix
Experts in Sound Science®

ECA
EUROPEAN COMPLIANCE ACADEMY

Prosonitron™ processing

Surface erosion compared with other ultrasonic devices

Device	Surface Intensity (W/cm²)	Process Intensity (W/litre)
Prosonitron P500 HD	0.1	~20
	1	~100
8mm Probe in 50ml	1	~20
	10	~100
Magnetostrictive 75mm Diameter Cell	10	~1000
	100	~10000


- Spreading of ultrasonic input over the whole surface of the cell allows significant specific power input to the fluid without demanding high surface intensity

© European Compliance Academy (ECA)

Prosonix
Experts in Sound Science®


Consideration of Ultrasonic Design Space

Understand zones of operation to match desired process duty



- Based on Prosonix field experience
- Residence time of several seconds required to allow time for process kinetics, heat transfer and mixing to keep pace
- Particular important consideration in crystallization

Sonocrystallization Design Space for QbD



Process understanding via:
Critical Quality Attributes and Critical Process Parameters

Prosonix Ultrasonic Particle Engineering 'Prosonitron™' Technology for Commercial Scale




SonoLab™ SL-250



5L Prosonitron™

© European Compliance Academy (ECA)



Prosonix
Experts in Sound Science®

Scale-Up Factors

Engineering Parameters to Maintain

- **Specific Power Input** – Intensity of acoustic power going into the process liquid (W / litre)
- **Specific Energy Input** – Intensity of acoustic energy used to achieve the process effect (J / litre)
- **Surface Intensity** – Acoustic power per unit area used to supply (W / cm²)
- **Residence Time** – Required duration inside field to achieve process effect (s)
- **Reynolds Number** – Internal mixing and flow pattern within acoustic field (-)

Prosonitron™ Design Benefits

- Ease of maintaining specific power input to the process liquid and surface intensity at all scales.
- Ability to reproduce flow regime and residence time by configuration choice.
- Result is a system that can deliver consistent processing throughout scale-up from pilot to production.

© European Compliance Academy (ECA)

Prosonix
Experts in Sound Science®



PAT & Complete Crystallization Control™

- Use CQA's and CPP's to understand and control manufacturing
- Use ultrasonic technology at scale to provide real time process control
- Need online measurement and feedback to achieve:
 - Control of ultrasonic energy
 - Control of temperature profiles
 - Turbidity / Lasentech for onset of nucleation
 - Particle size distribution
 - Polymorphic forms in solution
 - Supersaturation ?

© European Compliance Academy (ECA)



Prosonitron™ – Proven at commercial scale

"Best available crystallization technology" potential in pharma production

- Prosonitron™ at Pfizer, UCB & others...
- Pfizer Ireland acquire Prosonitron™ technology for Primary API Manufacture
- 4 year co-development relationship with UCB
- Many new & current trials at major and specialty pharma worldwide
- Primary for oral, but exciting developments in inhalation, nanosuspensions, parenteral, sub Q, and dermatological delivery



Prosonitron™ linked to Lasentec



- Prosonitron™ system in operation in alumina
- Over 3 years continuous service to date
- Follow on global deal with Alcoa

© European Compliance Academy (ECA)





Advanced Particle Engineering of Difficult to Manufacture products

e.g. Delivery by Inhalation

Product properties of particulate systems

Its much more than just size control....

Product property = f (Dispersity, Chemical Composition)

- Dispersity characterised by:
 - Particle Size
 - Aerodynamic \varnothing 0.5 - 5 μ m
 - Particle Shape
 - Conferred spherical geometry
 - Particle Surface Morphology
 - Minimise surface free energy - crystals
 - Particle Surface properties
 - Reduce contact area - spherical, rugged surface
- Control of interfacial interactions (adhesion/cohesion) is governed by surface forces
- Geometry, not surface chemistry, is the central design principle in controlling interfaces and their interactions.

Design of particles for end use properties!

Solution to Particle Technologies

Dr Ivan Marziano, Pfizer Material Science, APS Inhalation meeting 28 March 2007

Micronisation is often not feasible for inhaled medicines!


Solution to Particle technologies are based on the formation of highly supersaturated droplets either through spraying or dispersion in non-miscible media; Particle growth occurs within confined space (= droplet)

- Aerosol flow reactor (VTT)
- SAX™: solution atomisation and crystallization by sonication (**Prosonix**)
- Emulsion crystallization
- Quasi emulsion/spherical crystallization
- Cryogenic spray freezing/liquid extraction
- Spray freezing into liquids/spray freeze drying
- Spray drying
- EPAS: evaporative precipitation into aqueous solution (Dow)
- Segregated flow tubular reactor ("Bubbletube")

SAX™ – New "Solution to Particle" Technology


The future solution to ultimate particle control ?

- SAX™ builds on the Prosonitron™ IP
- Designed to produce 1 to 5 µm crystalline combination particles tailored uniquely for inhalation
- 20 customer studies completed to date
- Combination of proven unit operations:
 - Solution of API / mixtures of API + excipients
 - Atomization creates spherical droplets
 - Controlled evaporation
 - Controlled crystallization by Sonocrystallization
 - Simple and proven isolation procedure
- Tested on a range of APIs and NCEs
 - Stable crystalline particles
 - Combinations of 2 or 3 drugs possible
 - Improved *in-vitro* performance




ECA **Unprecedented Process & Product Control**
 EUROPEAN COMPLIANCE ACADEMY Optimal SAX™ drug particles should give better inhalation clinical performance

- Current destructive (i.e. micronisation) production techniques are severely limited
- API particles can be fully engineered from the the “ground up” via controlled SAX crystallization
- Improved SAX™ drug particles gives optimum formulation performance and patient benefit across all inhalation delivery platforms



Micronisation




SAX™

Size	<input checked="" type="checkbox"/>	Size	<input checked="" type="checkbox"/>
		Shape	<input checked="" type="checkbox"/>
		Surface Morphology	<input checked="" type="checkbox"/>
		Surface properties	<input checked="" type="checkbox"/>
		Crystallinity	<input checked="" type="checkbox"/>
		Stability	<input checked="" type="checkbox"/>

Prosonix
Experts in Sound Science®

© European Compliance Academy (ECA)

ECA **Seretide Combination: Single vs separate inhalers?**
 EUROPEAN COMPLIANCE ACADEMY Nelson et al J. Allergy Clin. Immunol. Vol 112 (1), 2003, 30



Raman laser analysis of Seretide metered-dose inhaler formulation on stage 4 ACI

Key:
 Fluticasone (green),
 Salmeterol (red)
 Co-association (yellow)


- 4 Separate clinical studies showed equivalence viz primary efficacy
- Consistent trend in favour of combination therapy
- Increased efficacy over concurrent use of same doses of same 2 drugs
- Co-deposition offers increased opportunity for synergistic interaction

Prosonix
Experts in Sound Science®

© European Compliance Academy (ECA)


ECA **SAX™ - Designed Combination Particles**
 EUROPEAN COMPLIANCE ACADEMY A new way to treat Asthma/COPD and other disease states?

- Eliminate variability associated with blending of 2 or more micronised powders
- Multiple API's in the correct dose ratio can be symbiotically crystallized in a single perfect particle
- Controlled Combination Therapy Delivery



00461 10KV X5,000 17mm

Fluticasone and Salmeterol (1:3.448)




10kV X7,000 2mm CEOS

Budesonide and Formoterol (1:17.71)

Prosonix
Experts in Sound Science®

© European Compliance Academy (ECA)

ECA **Improved Product Stability with SAX™**
 EUROPEAN COMPLIANCE ACADEMY Carrier based DPI formulation



Formulation	Time 0 (mcg)	40C/75RH (4 Weeks) (mcg)
Micronised FP	~25	~18*
Micronised SX	~2	~2*
SAX FP	~24	~26
SAX SX	~3	~3

- Carrier based DPI formulations (FP:SX 1:5)
- In vitro apparatus: NGI @ 60 L/min; Monohaler DPI device
- Seamless transition and increased stability

Prosonix
Experts in Sound Science®

© European Compliance Academy (ECA)

ECA Improved Dosing Compliance with SAX™
 Uniformity of FP-SX doses: lower stages of an NGI

Stage	Mic FP/SX (t=0)	SAX FP/SX (t=0)
Stage 2	~1.28	~0.98
Stage 3	~1.38	~0.95
Stage 4	~1.40	~0.93
Stage 5	~1.35	~0.92

The dose balance of micronised formulations of FP and SX not consistent between stages micronised batches.

Variability in the delivery of FP with respect to SX on the lower stages

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®

ECA Improved Dosing Compliance & Stability
 Uniformity of FP-SX upon storage

Stage	Mic FP/SX (4 wk, 40°C/75%RH)	SAX FP/SX (4 wk, 40°C/75%RH)
Stage 2	~1.38	~1.10
Stage 3	~1.18	~1.08
Stage 4	~1.08	~1.00
Stage 5	~0.98	~1.00

- Accelerated stability conditions affect the dose balance of micronised formulations
- SAX particles maintain high content uniformity of the actives

© European Compliance Academy (ECA) **Prosonix** Experts in Sound Science®



Summary

- Now possible to engineer and produce “Designed for Purpose” Particles
- QbD provides process understanding
- PAT provides repeatable and robust processes
- Key benefits include
 - Control particle size, shape, crystallinity, polymorphism
 - Improved batch consistency, filtration, isolation and drying
 - Replacement of problem physical seeding
 - Increased cGMP compliance
 - Improved formulation consistency, stability and performance
 - Increased return on investment
 - Reduced time to market
- Ultrasonic Particle Engineering and Sonocrystallization gives levels of control over isolation that current manufacture cannot deliver

© European Compliance Academy (ECA)



Contact

Christian Jones
BDM
Prosonix Ltd

+ 44 (0) 1865 784250

christian.jones@prosonix.co.uk

www.prosonix.co.uk

© European Compliance Academy (ECA)

